**SRRP**

**ENT**

**Septorhinoplasty - Septorhinoplastik**

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Dear Colleagues,

Have you also experienced the disappointment of searching for a particular instrument, only to end up wondering when you finally had it in your hands if it was the one you really wanted?

As a colleague in the field of ear, nose and throat surgery you will have been accustomed to a high standard of instruments in your clinical training and will be eager to maintain this standard in your own practice. Unfortunately, however, some of these standard instruments appear to have changed during the course of time. Even if an instrument is manufactured by the same company, it is frequently „not what it used to be“. Managing to find exactly the right instrument among the wide variety on offer is a tedious and time-consuming task which is aggravated by the fact that e.g. the „Cottle speculum“ you acquired is not the same as the one you had grown to appreciate in your many years of clinical practice.

Extra time and resources are again needed to have the instrument adapted to the desired and required expectations. In some cases subsequent modification is not even possible.

All the more welcome, then, is the news that Tekno Medical is doing something to counteract this trend. In cooperation with ENT surgeons, Tekno Medical is marketing instruments which have proved their reliability in many years of practice. Certain instruments are accompanied by a brief commentary outlining their outstanding features. In view of the minute dimensions involved in high-precision ENT work in body cavities, particular attention has been paid to the ultrafine design and slimline finish of the instruments.

In this catalogue series you will find individual sets of preselected instruments for the nose, paranasal sinuses and ear, which enable trouble-free performance of functional-aesthetic nose surgery, endonasal-endoscopic paranasal sinus surgery and middle-ear surgery including stapes surgery. The instruments have been selected on the basis of extensive operating experience in one of Germany's largest ENT clinics (Erlangen) and in my hospital practising.

If you have any questions or suggestions concerning the instruments, Tekno Medical and I are at your service at all times. Only a united effort will enable us to maintain and enhance our high standard of ENT surgery. Good instruments are an indispensable factor in achieving this aim. On this positive note, we wish you much enjoyment and pleasure working with these new instruments.

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
The complete set contains all listed items and a sterilization container (1/1, 500x280x150mm) with a wire basket.

- |    |              |   |    |              |  |
|----|--------------|---|----|--------------|--|
| ②③ | tk 56581-19W | Cottle Hammer, 19cm, 250gr, convex/flach<br><i>Cottle mallet, 19cm, 250gr, convex/flat</i>          | ③⑦ | tk 56645-20W | Saug-Raspatorium, 20cm, scharf<br><i>Suction raspatory, 20cm, sharp</i>                                  |
| ②④ | tk 56290-14W | Rubin Osteotom, 14mm, 18cm<br><i>Rubin Osteotome, 14mm, 18cm</i>                                    | ③⑧ | tk 56905-02W | Cottle Häkchen, 15cm, Fig.2<br><i>Cottle hooklet, 15cm, fig.2</i>  |
| ②⑤ | tk 56290-16W | Rubin Osteotom, 16mm, 18cm<br><i>Rubin Osteotome, 16mm, 18cm</i>                                    | ③⑨ | tk 56916-16W | Neivert Wundhaken, doppelendig, 20cm<br><i>Neivert retractor, double ended, 20cm</i>                     |
| ②⑥ | tk 56294-12W | Walter Osteotom, 12mm, 19cm<br><i>Walter Osteotome, 12mm, 19cm</i>                                  | ④① | tk 56919-13W | Aufricht Nasenhaken, 13,5cm<br><i>Aufricht nasal retractor, 13,5cm</i>                                   |
| ②⑦ | tk 56294-03W | Walter Osteotom, 3mm, 19cm<br><i>Walter Osteotome, 3mm, 19cm</i>                                    | ④② | tk 56750-00W | Cottle Daumenhaken<br><i>Cottle thumb hook</i>   |
| ②⑧ | tk 56292-03W | Walter Flachmeisel, Schwabenschwanz, 3mm, 19cm<br><i>Walter bone chisel, swallowtail, 3mm, 19cm</i> | ④③ | tk 56648-00W | Cottle Quetscher mit Klammer<br><i>Cottle bone crusher with clamp</i>                                    |
| ②⑨ | tk 56296-01W | Silver Osteotom, 18cm, gerade<br><i>Silver osteotome, 18cm, straight</i>                            | ④④ | tk 56755-00W | Winter Knorpelhalter<br><i>Winter cartilage holder</i>   |
| ③① | tk 56280-04W | Cottle Meißel, 4mm, 18cm<br><i>Cottle chisel, 4mm, 18cm</i>   | ④⑤ | tk 56185-18W | Bipolarepinzette, isoliert, 18cm, stumpf, baj.<br><i>Bipolar forcep, insulated, 18cm, blunt, bajonet</i> |
| ③② | tk 56669-21W | Aufricht Glabellaraspel<br><i>Aufricht Glabella rasp</i>  | ④⑥ | tk 56323-26W | Krause Nasenpolypenschnürer, 26cm<br><i>Krause nasal polypus snare, 26cm</i>                             |
| ③③ | tk 56751-15W | Diamond Nasenraspel, 19cm, gerade<br><i>Diamond nasal rasp, 19cm, straight</i>                      | ④⑦ | tk 56320-30W | Stahldrahtschlingen, 0,3mm, 12 Stück<br><i>Steel wire loops, 0,3mm, 12pcs.</i>                           |
| ③④ | tk 56686-03W | Fomon Raspel, TC, Fig. 3/4, 20cm<br><i>Fomon rasp, TC, fig. 3/4, 20cm</i>                           | ④⑧ | tk 56110-05W | Lösungsschale, 61x30mm, 0,063l, rostfrei<br><i>Bowl, 61x30mm, 0.063l, stainless</i>                      |
| ③⑤ | tk 56686-09W | Fomon Raspel, TC, Fig. 9/10, 20cm<br><i>Fomon rasp, TC, fig. 9/10, 20cm</i>                         |    | tk 56110-09W | Lösungsschale, 83x41mm, 0,16l, rostfrei<br><i>Bowl, 83x41mm, 0.16l, stainless</i>                        |
| ③⑥ | tk 56753-14W | Cottle Septummesser, 14cm<br><i>Cottle septum knife, 14cm</i>                                       |    | tk 56110-25W | Lösungsschale, 120x35mm, 0,25l, rostfrei<br><i>Bowl, 120x35mm, 0.25l, stainless</i>                      |
| ③⑦ | tk 56636-21W | Masing Periost-Elevatorium, 20cm<br><i>Masing periosteal elevator, 20cm</i>                         | ④⑧ | tk 56855-00W | Abziehstein<br><i>Oil stone</i>  |

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*Violent squeezing of the septum cartilage frequently causes destruction of the cartilage. This is often the reason why the inserted cartilage is no longer detectable after some time. In order to maintain cartilage consistency, working with a crushing forceps is recommended, since it largely saves the cartilage cells and ensures their survival.*

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*If the lateral nose walls are still too high after lateral and transversal osteotomy, further chiselling or rasping is no longer possible since the pyramid is too mobile. The only possibility left is cutting, and this bone-cutting forceps is extremely well suited for that procedure.*

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*Rasping the nasal prominence still propagated by many authors often causes traumatism of the nasal skin above the prominence and leads to skin atrophy. Instead, neat chiselling of unevennesses has proven effective. Afterwards, the sharp edges can be smoothed down very well using this diamond file.*

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*Placing an intranasal suture makes great demands on the surgeon's skill. This needle holder designed like a Blakesley permits fast working in the nose and is a wonderful help. The absence of a lock, allowing the opening and closing of the branches as the need arises is an additional advantage of this needle holder. Septum perforations can be closed very fast and nicely as well.*

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*Basically, the lifting hook has two functions, firstly, to keep open the cavity above the nasal dorsum during decollement and, secondly to protect the nasal skin from undesired perforation while ablating the nasal prominence. Many hooks available on the market are too short and narrow and, therefore, fail to meet the demands.*

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*The wire hook is a „must“ when working with open SRP without an assistant, which is often the case for attending doctors. This hook may be used nearly universally, and is very atraumatic due to its pointed hooks.*

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*This cartilage holder was specifically developed for working cartilage implants in nasal surgery. On the one hand, strongly deviated septa can be flattened and worked by means of scratching or sutures. On the other hand, sandwich implants can be sutured from several cartilage layers if a greater volume is required. It is also helpful for clamping thick costal cartilages and cutting away thinner slices of cartilage using an osteotomy blade.*

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*Sharp chisels identify the good nose surgeon. The chisels and osteotomes can be sharpened after use by means of this sharpening stone.*

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©2004 - TEKNO-MEDICAL Optik-Chirurgie GmbH & Co. KG
TEKNO is a registered trademark of Tekno-MEDICAL Optik-Chirurgie GmbH & Co. KG
**TEKNO-MEDICAL Optik-Chirurgie GmbH & Co. KG**