MOBI 2

MOBI - POWER TOOLS

MOBI - POWER TOOLS

PRÄZISION „MADE IN GERMANY“

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PRECISION „MADE IN GERMANY“

The new Power Tools series “MOBI 2” by Tekno-Medical guarantees a high performing motor system which will inspire you. Compared to the previous series “MOBI 1”, the new line is convincing through a more ergonomic design, an improved locking mechanism and a simplified preparation. All adapters and the PowerPack are compatible for the “MOBI 1” and “MOBI 2” system.

Unique durability of all components: Control electronics, motor and battery are in one removable PowerPack integrated. This means that only the motor housing of the respective power tool is autoclaved.

The high-quality processing of exclusive materials in combination with the user-friendly design assures comfortable and precise working. The intelligent positioning of the control elements allows an ergonomic one-hand operation.

Construction and design are made for the exigent everyday business of a hospital and ensure high reliability.

The TEKNO-MOBI system is your reliable partner in the field Of orthopedics and trauma surgery in the daily business of a hospital. The MOBI Power Tool is available in different configurations . All models have a quick connector included which ensures a simple and time saving replacement of the tools. The following models are available:

MOBI 2 DM Universal drilling / milling machine
MOBI 2 D Drilling machine
MOBI 2 O Oscillating saw
MOBI 2 RS Reciprocating / sternum saw
MOBI 2 ACL Reamer

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PowerPack

High performant LI-ION battery, without memory effect. The PowerPack may be used with all MOBI machines and recognizes automatically the Power Tool and adapts power, revolution and functions accordingly.

A six LED display indicates the battery capacity and function additionally as illumination during activation.

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 = Bohren  = Fräsen

Drilling and milling function (only MOBI DM)

By turning the ring, the drilling or milling function may be selected.

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MOBI D / DM – Drilling/milling machines

All units main functions may be controlled by three elements. A sliding switch is used to set the sense of rotation. With the sliding switch in middle position, the Power Tool may not be activated. The center knob regulates power/speed of the Power Tool. A switch allows the selection of permanent or oscillating mode.

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LONG-LIVING CONSTRUCTION

The PowerPack is in fact, the center-piece of the MOBI Power Tools. It is inserted into the handle of the machine. Beside the power-battery this module also consists of the impulse- as well as the processor unit. Thanks to the latest LI-ION technology a high capacity battery without memory-effects is provided. The PowerPack can be inserted in all machines available. The integrated processor recognizes automatically in which Power Tool it is used and adapts its power, revolution and functions according to the appropriate machine/insert.

Treatments by the autoclave are a high strain for batteries and electric components what can have negative effects on the life expectancy. This fact has been considered in the construction of the Tekno-Medical MOBI system. For this reason all sensitive components have been included into the TEKNO PowerPack, which is not autoclavable.

Due to this two-piece structure and high-quality manufacturing, the MOBI Power Tools achieve a longevity as well as a constant efficiency over years.

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MOBI D (tk 38501-00)

- Dimensions hand piece (without drill bit):
140 x 111 x 206 mm
- Weight hand piece with PowerPack: 1660 g
- Infinitely variable speed:
0 - 1000 rpm
- Cannulation: Ø 4.3mm
- Protecting class: B, EN 60601-1

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MOBI DM (tk 38501-01)

- Dimensions hand piece (without drill bit):
166 x 111 x 207 mm
- Weight hand piece with PowerPack: 1850 g
- Infinitely variable speed:
0 - 1000 rpm (Mode Drill)
0 - 250 rpm (Mode Ream)
- Cannulation: Ø 4.3mm
- Protecting class: B, EN 60601-1

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MOBI ACL (tk 38501-09)

- Dimensions hand piece (without saw blade):
140 x 111 x 206 mm
- Weight hand piece with PowerPack: 1550 g
- Infinitely variable speed:
0 - 9000 rpm (Step I)
0 - 11000 rpm (Step II)
- Protecting class: B, EN 60601-1

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MOBI O (tk 38501-02)

- Dimensions hand piece (without saw blade):
166 x 111 x 211 mm
- Weight hand piece with PowerPack: 1760 g
- Infinitely variable speed:
0 - 9000 rpm (Step I)
0 - 11000 rpm (Step II)
- Protecting class: B, EN 60601-1

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MOBI RS (tk 38501-04)

- Dimensions hand piece (without saw blade):
185 x 111 x 206 mm
- Weight hand piece with PowerPack: 1750 g
- Infinitely variable speed:
0 - 7500 rpm (Step I)
0 - 10000 rpm (Step II)
- Protecting class: B, EN 60601-1

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PowerPack (tk 38501-05)

- Dimensions: 89 x 87 x 102 mm
- Weight: 760 g
- Type: LI-ION
- Max. voltage: 16.8 V
- Nominal voltage: 14.4 V
- Capacity: 2.1 Ah
- Typical charging time: < 90 min.

tk 38501-18

AO-Aufsatz, klein
AO chuck, small

tk 38501-16

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tk 38501-26

DIN-Aufsatz
DIN chuck

tk 38501-21

Jacobs Bohrfutter, klein (bis 4 mm Ø)
Jacobs drill chuck, small (up to 4 mm Ø)

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
MOBI - DRILL / MILLING MACHINE

tk 38501-11

Spickdrahtfutter, 1-4 mm
Wire/pin driver, 1-4 mm

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