HF-GERÄT TOM 500

HF-UNIT TOM 500

TEKNO TOM 500

TEKNO TOM 500

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Along with the steadily growing requirements by OR and offices comes an increasing growth of application programs, the HF units and their handling get more complicated.

Tekno-Medical has set its goal for the new TOM 500, to portray all important data and buttons as simple as possible for a quick and intuitive handling of the units at any time.

All relevant information at that time and all command elements are shown on the 7" (17,8cm) touch-display. The clear depiction is the base for safe reading and handling and the avoidance of stress situations.

The conception of the Tekno TOM 500 units impresses by their universal possibilities for use. The units are reliable partners for all interventions in the open surgery as well as in the endoscopic surgery. Tekno TOM 500 units are optionally available with an integrated Argon-module. In combination with the smoke evacuators Tekno-Medical can offer a complete system.

Experience of many years and first class engineering of Tekno-Medical build the base for the development of the TOM 500 units, which stand out due to their performance, security technology and the excellent handling.

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OPERATION


Easy menu-driven operation through three levels

The user reaches the program levels or the configuration levels easily from the base level.

The program levels, where the user can see settings of the chosen program, can be reached by touching the individual fields.

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
PROGRAM LEVEL

Easy setting of programs and performance.

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
STATUS LEVEL

On the large 7" touch-display the user has all relevant information in the point of view. The Tekno TOM 500 units show besides the chosen programs, the set performance and the selected neutral electrode including the contact to the patient.

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CONFIGURATION LEVEL

The user can save various presetting's, to customize the HF-unit.

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PROGRAMMS

The user has 99 free definable programs available. These programs are already set for different applications.

Personalizing the programs, e.g. by an individual assignable name, is self-evident for the Tekno TOM 500 units.


The following are available as standard modes:

- SEALING / AUTO SEALING function – for sealing vessels up to Ø 8mm
- Bipolar TUR-P / TUR-B / ARTHRO CUT MODI
- 7 different modes (5 monopolar / 2 bipolar) with 10 effects in CUT-mode each
- Specific modes for flexible endoscopy
AUTO PAPILO PURE CUT
AUTO POLIPO BLEND CUT
AUTO ENDO CUT
- PURE CUT PULSED / BLEND CUT PULSED for open surgery and flexible endoscopy
- Auto start/stop function
- Simultaneous working with 2 instruments (monopolar field)

Each unit is provided with an automated program for the control and adjustment of the output currents with a PER (Power Efficiency Rating) of more than 98%. The system automatically adopts all currents according to the specific tissue features with two different systems according to the selected mode:

ADC System – Self adjusting with constant power

APC System – Self adjusting with constant voltage

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SAFETY

The system also compares the output power with the set values, to In the stage of developing, Tekno-Medical has paid particular attention to the safety of the HF-units, next to user-friendly use application performance. All relevant functions for safety, e.g. regulation of the output power control of the neutral electrode, etc. are permanently monitored in the background by microprocessors, without disturbing the user.

By using REM neutral electrodes (Return Electrode Monitoring), the unit permanently checks the contact between patient and neutral electrode. Is the contact area reducing, the unit is giving an acoustic and optical alarm and decreases the output power automatically step by step to zero, to prevent burns on the patient. When using small neutral electrodes, the unit limits the maximum output power automatically (baby function).

Beside the standard neutral plate connection, the Tekno TOM 500 disposes additionally of a connection for neutral plates with Valleylab™ plugs.

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MODEL VARIATIONS

Tekno TOM 500 units are available in 4 different variations. The user has the opportunity to choose an HF-unit for a good value for money which adopts perfectly to his needs.

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The footswitch is equipped with an additional selector button, allowing the user to choose between different pedal selections. The TOM 500 units may as well be commanded the classical way by two separate foot-switches.

Universal (Mono-/Bipolar)		Bipolar
tk 90016-01	tk 90016-02	tk 90017-00
Standard Ausführung <i>Standard version</i>	Breite Ausführung <i>Wide version</i>	Standard Ausführung <i>Standard version</i>

	tk 90295-01	tk 90295-02	tk 90302-11
<i>for units</i>	1	2	1 or 2
<i>Brakes</i>	2	2	2
<i>Antistatic wheels</i>	x	x	x
<i>Holder for gas bottles</i>	x (optional)	x (optional)	x

ZUBEHÖR

ACCESSORIES

STANDARD SETS

For TEKNO TOM 500		Standard		Argon		
		tk 90010-51	tk 90010-61	Basic	LAS	Gastro
				tk 90011-53	tk 90011-54	tk 90011-55
Handle insertion in mm		Ø 2.4	Ø 4.0			
Double pedal switch, UNIVERSAL	tk 90016-01	x	x	x	x	x
Electrode Handle with manual control	tk 90020-00	x		x	x	
Electrode Handle with manual control	tk 90030-01		x			
Argon Electrode Handle with manual control	tk 90301-10			x	x	
Electrode Handle without manual control	tk 90021-00	x		x	x	
Cable for flexible electrodes	tk 90301-18					x
Knife electrodes (2x)	tk 90022-01	x		x	x	
Knife electrodes (2x)	tk 90032-01		x			
Knife electrode, rigid, 40 mm	tk 90301-16			x	x	
Knife electrodes, 100mm	tk 90022-05	x		x	x	
Rhombic knife electrode	tk 90032-02		x			
Needle electrode, thick	tk 90023-01	x		x	x	
Needle electrode, thick	tk 90033-01		x			
Needle electrode, thin	tk 90023-02	x		x	x	
Needle electrode, thin	tk 90033-03		x			
Ball electrode, Ø 2.5 mm	tk 90024-01	x		x	x	
Ball electrode, Ø 4.0 mm	tk 90034-01		x			
Ball electrode, Ø 4.0 mm	tk 90024-02	x		x	x	
Ball electrode, Ø 6.0 mm	tk 90034-02		x			
Loop electrode, Ø 5.0 mm	tk 90025-01	x		x	x	
Loop electrode, Ø 10.0 mm	tk 90035-01		x			
Coagulation electrode, rigid, 25 mm	tk 90301-11			x	x	
Coagulation electrode, rigid, 100 mm	tk 90301-12			x		
Coagulation electrode, rigid, 320 mm	tk 90301-13				x	
Flexible electrode, Ø 2.3 mm, 1000 mm	tk 90301-20					x
Disposable neutral electrode with cable, sterile	tk 90029-16	x	x	x	x	x
Tube for gas input with quick-action coupling	tk 90302-14			x	x	x
Pressure reducer for argon gas bottles	tk 90302-13			x	x	x
Filter for argon gas output	tk 90302-24			x	x	x
Luerlock adaptor for argon gas filter	tk 90302-25			x	x	x

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© 2015 – TEKNO-MEDICAL Optik-Chirurgie GmbH
TEKNO is a registered trademark of TEKNO-MEDICAL Optik-Chirurgie GmbH