

Saug- / Spülpumpe MULTI PUMP

Suction- / irrigation unit MULTI PUMP



MULTI PUMP

tk 750-5060

Technische Daten

- Druckbereich 0 - 750 mmHg
- Spüleistung 0 - 3,0 l/min
- Saugleistung 0 - 3,0 l/min

- Stromversorgung 100 / 240 VAC
- Frequenz 50 - 60 Hz
- Max. Leistungsaufnahme 36 W
- Gewicht ca. 7 kg
- Abmessung (B x T x H) 304 x 345 x 115 mm

Im Lieferumfang enthalten

- Silikonschlauchsatz, komplett, unsteril, autoklavierbar (tk 750-5061)
- Kanüle, 5 mm, autoklavierbar (tk 750-5065)
- Kanüle, 10 mm, autoklavierbar (tk 750-5066)
- Adapterrohr, unsteril, autoklavierbar (tk 750-5064)
- Handsteuerung pneumatisch (tk 750-5062)
- Steuerungsschlauch, autoklavierbar (tk 750-5063)

Technical data

- Pressure range 0 - 750 mmHg
- Irrigation rate 0 - 3.0 l/min
- Suction rate 0 - 3.0 l/min

- Power supply 100 / 240 VAC
- Frequency 50 - 60 Hz
- Max. power consumption 36 W
- Weight approx. 7 kg
- Dimensions (W x D x H) 304 x 345 x 115 mm

Included in delivery

- Silicone tube set, complete, not sterile, autoclavable (tk 750-5061)
- Canula, 5 mm, autoclavable (tk 750-5065)
- Canula, 10 mm, autoclavable (tk 750-5066)
- Adaptor tube, not sterile, autoclavable (tk 750-5064)
- Hand control pneumatic (tk 750-5062)
- Control tube, autoclavable (tk 750-5063)