

## TEKNO LIZARD TEKNO LIZARD



Ø 10.0 mm - 330 mm

**tk 707-090**

**tk 707-091**

Ersatzmesser (Pack á 5 Stück)  
 Replacement knife (Box of 5 pieces)

Das Instrument für die zuverlässige und schnelle Gefäßversiegelung. Die Bauweise als Kombi-Instrument erlaubt das sichere Versiegeln von Blutgefäßen und Gewebe sowie anschließenden, mittigen Schnitt. Dies erspart dem Operateur den aufwändigen Wechsel der Instrumente sowie die Positionierung für den Schnitt. Die Wiederverwendbarkeit des Instruments erlaubt einen kostengünstigen Einsatz. Die Klinge wird nur einmal verwendet und lässt sich vor der Aufbereitung einfach und schnell einsetzen um gleichbleibende optimale Ergebnisse zu garantieren.

*The instrument for fast and secure vessel sealing. The design as an instrument combination allows the secure sealing of blood vessels and tissue followed by a centered cut. Thus the surgeon does not need to change the instruments and search again for the correct cutting position. The cost-efficient instrument is reusable. The blade is a single-use component and is easy to fix before the intervention to grant a perfect cut each time.*