


Flaches Laryngoskop, HD, lateral beleuchtet

Flat Laryngoscope, HD, lateral illumination

autoklavierbar / autoclavable
 134°C / 273°F - 18 min.



Reizfrei durch extra flaches Design
 Less irritating due to flat design



Ø d1 = 7,2mm, d2 = 12,4 mm, 183,0 mm	tk 700-047*	tk 700-046*
--------------------------------------	-------------	-------------

- Ø, Arbeitslänge
- Ø, Working length
- Ø, Longueur de travail
- Ø, Longitud de trabajo
- Ø, Lunghezza operativa

