



ZERLEGBARE ARTHROSKOPIE-STANZEN  

---

*DISMANTABLE ARTHROSCOPY PUNCHES*

## ZERLEGBARE ARTHROSKOPIE-STANZEN


### DISMANTABLE ARTHROSCOPY PUNCHES


Durch Knopfdruck lässt sich der Schaft öffnen und ermöglicht so eine optimale Reinigung und Aufbereitung der Stanze. Dank der Drehbarkeit des Schaftes von fast 360° lässt sich die Stanze problemlos im Siebkorb positionieren. Bei geöffnetem Schaft sind beide Teile noch über einen Stift miteinander verbunden, wodurch verhindert wird, dass Einzelteile während des Aufbereitungsprozesses verwechselt werden oder verloren gehen. Intuitive und einfachste Handhabung runden das Produkt ab.

At the push of a button, the sheath will be open. This enables the possibility for an optimal cleaning and processing. Thanks to the rotation of almost 360° of the sheath, the punch can be placed smoothly into the wire basket. In opened condition, both parts are still connected to each other, which avoids loss and mistakes during the cleaning and sterilization process. Intuitive and simple handling complete this product.




#### DUCKBILL MEDIUM


	x = 2.5 mm	y = 3.2 mm	<b>tk 39400-10Z</b>
---	------------	------------	---------------------

	x = 2.5 mm	y = 3.2 mm	<b>tk 39400-14Z</b>
---	------------	------------	---------------------

Korbstanze  
*Basket punch* gerade  
*straight*




#### OVAL SHAPE


	x = 1.7 mm	y = 4.4 mm	<b>tk 39400-60Z</b>
---	------------	------------	---------------------

Korbstanze  
*Basket punch* gerade  
*straight*





#### LOW PROFILE

	x = 1.9 mm	y = 3.2 mm	<b>tk 39400-20Z</b>
---	------------	------------	---------------------

	x = 1.9 mm	y = 3.2 mm	<b>tk 39400-26Z</b>
---	------------	------------	---------------------

Korbstanze  
*Basket punch* gerade  
*straight*



	x = 2.0 mm	y = 1.6 mm	<b>tk 39660-02Z</b>	<b>tk 39661-02Z</b>
---	------------	------------	---------------------	---------------------

Korbstanze  
*Basket punch* rechts  
*right* links  
*left*



#### ALLIGATOR

a = 3.4 mm	<b>tk 39400-95Z</b>
------------	---------------------

Fasszange  
*Grasping forceps* gerade  
*straight*

© 2018 – TEKNO-MEDICAL Optik-Chirurgie GmbH  
TEKNO is a registered trademark of TEKNO-MEDICAL Optik-Chirurgie GmbH  
Printed in Germany 11/2018 – Nachdruck, auch auszugsweise, verboten