




SIEBKORB FÜR RESEKTOSKOPIE-SET

BASKET FOR RESECTOSCOPY SET

SIEBKORB FÜR RESEKTOSKOPIE SET

BASKET FOR RESECTOSCOPY SET

ALL IN ONE

Durch die speziell platzierten Silikonhalterungen, findet ein komplettes Resektoskopie Set in nur einem Siebkorb Platz. Dies sorgt durch eine übersichtliche und schnellere Zuordnung für Zeitersparnis im Klinik-Alltag.


Dank der flexiblen und universellen Halterungen kann der Siebkorb sowohl für monopolare, als auch bipolare Sets in verschiedenen Größen (22/19 Charr. und 26/24 Charr.) verwendet werden.

Zu verwenden ist der Siebkorb mit einem 1/1 Container:

Höhe 100mm für 1 Siebkorb

Höhe 200mm für 2 Siebkörbe

Höhe 260mm für 3 Siebkörbe



540 x 245 x 75 mm

tk 85109-05

Lochblechsiebkorb mit Deckel
Perforated plate basket with lid



CLASSIC **1/1**

No

FILTER

2 x

FILTER



585 x 274 x 100 mm

tk 87170-B0F-SI

tk 87170-B1F-SI

x 200 mm

tk 87200-B0F-SI

tk 87200-B1F-SI

x 260 mm

tk 87210-B0F-SI

tk 87210-B1F-SI

ALL IN ONE

Due to the specially placed silicone holders, a complete resectoscopy set can be stored in just one basket.

This saves time in clinical workaday life due to a clear and faster allocation.

Thanks to the flexible and universal holders, the basket can be used for monopolar and bipolar sets in different sizes (22/19 Charr. and 26/24 Charr.).

The basket is suitable for 1/1 containers:

Height 100mm for 1 basket

Height 200mm for 2 baskets

Height 260mm for 3 baskets




Abbildung zeigt Beispiel mit Resektoskopie-Set 24/26 Charr., passiv, bipolar. Lieferumfang ohne Instrumente.

Image shows example for resectoscopy set 24/26 Charr., passive, bipolar. Scope of delivery without instruments.



2 x

FILTER



585 x 278 mm

tk 87170-L1F-SI

585 x 278

tk 87170-L1F-YE

585 x 278

tk 87170-L1F-GR

585 x 278

tk 87170-L1F-BL

585 x 278

tk 87170-L1F-RE

585 x 278

tk 87170-L1F-BK

© 2020 – TEKNO-MEDICAL Optik-Chirurgie GmbH

TEKNO is a registered trademark of TEKNO-MEDICAL Optik-Chirurgie GmbH

Printed in Germany 10/2020 – Nachdruck, auch auszugsweise, verboten