



LEICHTGEWICHT NADELHALTER

LIGHTWEIGHT NEEDLE HOLDER

LEICHTGEWICHT NADELHALTER

LIGHTWEIGHT NEEDLE HOLDER

Die Tekno Leichtgewicht Nadelhalter mit geradem oder gebogenem Axialgriff überzeugen durch ihr ergonomisches Design, ihre Funktionalität (Fassen verschiedenster Nadelgrößen) sowie ihre Langlebigkeit Dank geringer Abnutzung der Maulteile.

Aufgrund des geringen Gewichtes (Handgriffe aus Aluminium) wird ein langes und ermüdungsfreies Arbeiten ermöglicht.

Des Weiteren besitzen alle Leichtgewicht Nadelhalter einen integrierten Überlastschutz, welcher das Produkt vor Beschädigungen durch zu hohe Krafteinwirkung schützt.

Alle Modelle sind mit einem LuerLock Spülanschluss ausgestattet, welcher eine optimale Reinigung und Wiederaufbereitung garantiert.

Zur Auswahl stehen diverse Modelle in verschiedenen Längen:


- Ø 5mm als nicht zerlegbares Komplett-Instrument oder als zerlegbare Variante (Griff, Schaftrohr und Einsatz).
- Ø 3mm als nicht zerlegbares Komplett-Instrument.

The light weight needle holders with straight or curved Axial handle convincing through its ergonomic design, the function (grasp for different needle sizes), thanks to the low wear of the jaw parts, with a long resistance.

The light weight (handles made of aluminium) enables a long and fatigue-proof working process. Furthermore all light weight needle holders have an integrated overload safety system, which protects the product from damage made through force effect. All models are equipped with a Luer-Lock flushing port, which guarantees for an optimal cleaning and reprocessing.


There are different models in different lengths to choose from:

- Ø 5mm as not detachable complete instrument or as detachable version (handle, shaft and insert).
- Ø 3mm as not detachable complete instrument.



LEICHTGEWICHT NADELHALTER, NICHT ZERLEGBAR

LIGHTWEIGHT NEEDLE HOLDER, NON DETACHABLE




		Ø 3.0 mm	Ø 3.0 mm	
	Nadelhalter, TC, Maul gerade <i>Needle holder, TC, jaw straight</i>	250 mm 330 mm	tk 708-560 tk 708-562	tk 708-561 tk 708-563
	Nadelhalter, TC, Maul links gebogen <i>Needle holder, TC, jaw left curved</i>	250 mm 330 mm	tk 708-564 tk 708-566	tk 708-565 tk 708-567
	Nadelhalter, TC, Maul rechts gebogen <i>Needle holder, TC, jaw right curved</i>	250 mm 330 mm	tk 708-568 tk 708-570	tk 708-569 tk 708-571

		Ø 5.0 mm	Ø 5.0 mm	
	Nadelhalter, TC, Maul gerade <i>Needle holder, TC, jaw straight</i>	330 mm 450 mm	tk 708-540 tk 708-542	tk 708-541 tk 708-543
	Nadelhalter, TC, Maul links gebogen <i>Needle holder, TC, jaw left curved</i>	330 mm 450 mm	tk 708-544 tk 708-546	tk 708-545 tk 708-547
	Nadelhalter, TC, Maul rechts gebogen <i>Needle holder, TC, jaw right curved</i>	330 mm 450 mm	tk 708-548 tk 708-550	tk 708-549 tk 708-551
	Nadelhalter Flamingo, diamantiert, Maul nach oben gebogen <i>Needle holder Flamingo, diamond coated, jaw curved up</i>	330 mm 450 mm	tk 708-552 tk 708-554	tk 708-553 tk 708-555
	Nadelhalter 1x2 Zähne, TC, Maul gerade <i>Needle holder 1x2 teeth, TC, jaw straight</i>	330 mm 450 mm	tk 708-556 tk 708-558	tk 708-557 tk 708-559

Axial Handgriff
gerade
*Axial handle
straight*

Axial Handgriff
gebogen
*Axial handle
curved*

LEICHTGEWICHT NADELHALTER, ZERLEGBAR
LIGHTWEIGHT NEEDLE HOLDER, DETACHABLE



Nadelhalter, TC, Maul gerade
Needle holder, TC, jaw straight

330 mm

Ø 5.0 mm
tk 708-640

Ø 5.0 mm
tk 708-641

Ø 5.0 mm
tk 708-634



Nadelhalter, TC, Maul links gebogen
Needle holder, TC, jaw left curved

330 mm

tk 708-644

tk 708-645

tk 708-635



Nadelhalter, TC, Maul rechts gebogen
Needle holder, TC, jaw right curved

330 mm

tk 708-648

tk 708-649

tk 708-636



Nadelhalter Flamingo, diamantiert,
 Maul nach oben gebogen
*Needle holder Flamingo, diamond coated,
 jaw curved up*

330 mm

tk 708-652

tk 708-653

tk 708-637



Nadelhalter 1x2 Zähne, TC, Maul gerade
Needle holder 1x2 teeth, TC, jaw straight

330 mm

tk 708-656

tk 708-657

tk 708-638



Nadelhalter, Maul gerade, gerieft
Needle holder, jaw straight, serrated

330 mm

tk 708-660

tk 708-661

tk 708-639

Komplett-
 instrument
*Complete
 instrument*

Komplett-
 instrument
*Complete
 instrument*

Einsatz allein
 Insert only



tk 708-630

Axial Handgriff, gerade
Axial handle, straight



tk 708-631

Axial Handgriff, gebogen
Axial handle, curved



330 mm

Ø 5.0 mm

tk 708-633

Metallschaft
Metal shaft