


ISO-FOCAL ZOOM ADAPTER FÜR 4K KAMERA-SYSTEME

ISO-FOCAL ZOOM ADAPTER FOR 4K CAMERA SYSTEMS



f = 16-32 mm

tk 754-3914 UHD

ISO-FOCAL ZOOM Objektiv
ISO-FOCAL ZOOM Objective

ISO-Focal Zoom Objektiv mit brillianter Bildqualität in Verbindung mit 4K Kamera-Systemen. Der spezielle TV-Adapter bleibt fokussiert während Sie den Zoom und damit die Brennweite ändern. Somit entfällt in der Praxis das lästige Nachfokussieren des Bildes. Ideal für alle endoskopischen Bereiche. Der TV-Adapter verfügt über einen C-Mount Anschluss sowie einem Fingerrad zur Brennweitenregulierung und einem Ring zur Fokussierung. Der Adapter ist einlegbar.

ISO-focal zoom objective with brilliant image quality in combination with 4K camera systems. The special TV-adapter with zoom keeps focussed while changing the focal length. The major advantage is that there will be no more time lost due to re-focussing the image after changing the zoom. Perfect for all endoscopic fields.

The TV-adapter has a C-mount connection, a double control ring for focal length and focus. The adapter is soakable.

© 2020 – TEKNO-MEDICAL Optik-Chirurgie GmbH
TEKNO is a registered trademark of TEKNO-MEDICAL Optik-Chirurgie GmbH
Printed in Germany 10/2020 – Nachdruck, auch auszugsweise, verboten / 0118